




## MAIN FEATURES

- > Adjusted for use in cytology, urology, microbiology, virology, oncology
- > Wide range of consumables for different applications
- > Centrifugation of up to 12 samples
- > Robust design
- > Microprocessor controlled functions of cover, time and speed
- > Adjustable acceleration and braking system
- > Special START - STOP features
- > FAST SPIN function at models CytoCENTRIC 150, CytoCENTRIC 250 and CytoCENTRIC 350
- > Possibility of programming



CYTO CENTRIFUGE	CytoCENTRIC 150 Code 455.000/TH4 Code 455.000/TH8	CytoCENTRIC 322 Code 464.000/TH	CytoCENTRIC 250 Code 472.000/TH	CytoCENTRIC 350 Code 473.000/TH
No. of places in rotor	4 or 8	12	8	12
Max. speed (RPM)	2.000	2.000	2.000	2.000
Brushless motor	NO	NO	YES	YES
Noise level (dBA)	54	56	48	49

## ACCESSORIES FOR CytoCENTRIC CENTRIFUGES

ROTOCLIPS STAINLESS STEEL JC 390 Part No. 312491	ROTOFUNELL CHAMBERS JC 406 Part No. 312488	FILTER CARDS JC 407 Part No. 309097	SLIDES JC 311 Part No. 303369	ECOfunnel® Reusable JC 313 Part No. 334724	ECOseal JC 315 Part No. 334729	ECOslides JC 314 Part No. 334733
4 pcs in a pack	4 pcs in a pack	200 pcs in a pack	100 pcs in a pack	4 pcs in a pack	100 pcs in a pack	100 pcs in a pack
	 <p>For sample volumes of up to 0.5ml Sedimentation field with a diameter of 6mm</p>			 <p>For preparing large-volume samples of up to 6ml Sedimentation field of 22x15mm for a higher cell yield Only for rotor with 12 places</p>		



We reserve the right to alter specification details etc. without prior notice or liability!

# DOMEL®

DOMEL, d.o.o., Otoki 21, 4228 Železniki  
 BU Laboratory Systems, Na Plavžu 79, 4228 Železniki  
 Phone: +386 4 5117 500 Fax: +386 4 5117 501  
 WEB: <http://www.domel.com/products/laboratory-equipment>

E-mail: [info@tehnica.si](mailto:info@tehnica.si)